

ANIMAS PUMP UPLOADING FOR PATIENTS

Patients need the following to upload the Animas pump at home:

1. Internet access, mac or PC
2. Diasend account, which is set up at www.diasend.com
3. Device cable with infrared dongle (can be ordered through 1-888-DIASEND or www.animas.com/get_software1), blue landing pad is optional
4. Diasend Uploader software (patients are directed to install when registering)

To register:

1. Go to www.diasend.com/register
2. Establish a username (typically the patient's email) and password
3. Device serial number (on back of pump)
4. Clinic id: 94-39631 (TDE)
5. Install Uploader software (follow prompts)

Uploading directions once patient is registered and the software is installed:

1. Double click on Diasend Uploader icon on PC desktop or through applications folder on a mac
2. A window will appear once Uploader is opened with the message "scanning for device"
3. Connect the IR dongle with the USB plug end to the computer
4. Suspend pump (the pump screen needs to be "awake"/illuminated for upload to occur)
5. Face pump infrared sensor window (back of pump, left of serial number) toward center of IR on dongle. The blue landing pad can be used but isn't essential.
6. The uploading process will usually take 5+ minutes
7. Once done, patients can select among 3 options: "quit and go to diasend.com" (to view reports) , "upload another device" or "quit".

Demo videos on uploading for patients available at www.diasend.com/us/patient

FAQ at www.diasend.com/us/support/faq-contact

