




## Registration and Upload Instructions

- 1) Go to [www.Carelink.Minimed.com](http://www.Carelink.Minimed.com)
- 2) Click on the **Sign Up Now** button
  - Choose your language
  - Accept the “Terms of Use and Privacy Statement”
  - Complete the Enrollment form.
    - i. **Username:** your first initial and last name.
      1. Example: alincoln
      2. Or full name and last name if that username is taken. Example: abrahamlincoln
    - ii. **Password:** dr\_\_\_\_ Example: drseuss



- 3) Return to the Carelink Personal main page and sign in using your newly created username and password.
- 4) Once signed in, choose “Upload data from my Device.”
  - Allow permission for CareLink to install the software “Java” if you are prompted
  - Choose the device you want to upload:



- 5) You will be asked to **Check Pump Status** and then **Identify the pump...** with your pump serial number.

- 6) You will then Choose your link Device
  - This is when you will choose your Contour Next Link meter

**Note:** Wait until prompted to Connect the meter to your computer

- Select **Next >**
- 7) You will now connect your meter to your computer and select **Finish**
- 8) Stay close to your computer while your pump is uploading. Once completed, you may log out or review your reports.



### HAVING TROUBLE?

Visit [www.medtronicdiabetes.com/carelinksetup](http://www.medtronicdiabetes.com/carelinksetup) for step-by-step instructions or call the Medtronic 24-Hour HelpLine anytime at 1.800.646.4633, option 1 for assistance.