

T-SLIM PUMP UPLOADING FOR PATIENTS

To upload the T-Slim pump patients need the following:

- USB Cable, pump, computer w/ internet access, memory card (optional)
- Mac or PC
- t:connect software

The following link will take the patients through the process of uploading the software and creating an account:

<https://tconnect.tandemdiabetes.com/Login.aspx?ReturnUrl=%2f>

Once patients have installed the software and created an account, they follow the next steps to upload:

1. Sign in at <http://tconnect.tandemdiabetes.com> w/ username and password
2. Connect pump to computer using USB cable
3. Open t:connect software
4. Select “upload “
5. Once they see the “upload successful” message, patients can either select “view reports” or “save reports”

